

Anti-MNX1/HLXB9 Antibody (5J903)

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse)
Molecular Weight:	Theoretical: 41 kDa. Actual: 48 kDa.
Clone:	5J903
Purification:	Protein A purified

Applications

Verified Activity:	25 µg total protein per lane of various lysates (see on figure) probed with MNX1/HLXB9 monoclonal antibody, unconjugated (TMAB-08893) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IP
Recommended	WB: 1:500-2000; IP: 1:20-50

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human MNX1
Antigen Species:	Human
Gene ID:	3110
Uniprot ID:	P50219

Research Background

Putative transcription factor involved in pancreas development and function.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481