

## Anti-MSL3L2 Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 51 kDa. Actual: 51 kDa.
Purification:	Protein A purified

### Applications

Verified Activity:	Sample: Hcclm3 (Human) Cell Lysate at 30 µg Primary: Anti-MSL3L2 (TMAB-09039) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 51 kD Observed band size: 51 kD
Application:	WB
Recommended	WB: 1:500-2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human MSL3L2
Antigen Species:	Human
Gene ID:	151507
Uniprot ID:	P0C860

### Research Background

MSL3L2 contains 1 chromo domain. It may be involved in chromatin remodeling and transcriptional regulation.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481