

## Anti-MST1 Antibody (5L825)

### Product Details

|                   |                                      |
|-------------------|--------------------------------------|
| Ig Type:          | IgG                                  |
| Reactivity:       | Human,Mouse,Rat                      |
| Molecular Weight: | Theoretical: 80 kDa. Actual: 80 kDa. |
| Clone:            | 5L825                                |
| Purification:     | Protein A purified                   |

### Applications

|                    |  |
|--------------------|--|
| Verified Activity: | 25 µg total protein per lane of various lysates (see on figure) probed with MST1 monoclonal antibody, unconjugated (TMAB-09042) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min. |
| Application:       | WB   |
| Recommended        | WB=1:500-2000  |

### Properties

|                      |   |
|----------------------|---|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping:            | Shipping with blue ice.   |

### Antigen Details

|                  |                                   |
|------------------|-----------------------------------|
| Immunogen:       | A synthesized peptide: human MST1 |
| Antigen Species: | Human                             |
| Gene ID:         | 4485                              |
| Uniprot ID:      | P26927                            |

### Research Background

Probably has no proteolytic activity, since crucial characteristic of serine proteases catalytic sites are not conserved.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481