

Anti-MUC16 Antibody (1E460)

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 1519 kDa. Actual: 180-350 kDa.
Clone:	1E460
Purification:	Protein A purified

Applications

Verified Activity:	25 ug total protein per lane of various lysates (see on figure) probed with MUC16 monoclonal antibody, unconjugated (TMAB-09094) at 1: 1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,ICC/IF,FCM
Recommended	WB: 1:500-2000; ICC/IF: 1:50-200; FCM: 1:50-100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Mucin 16
Antigen Species:	Human
Gene ID:	94025
Uniprot ID:	Q8WXI7

Research Background

Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481