

Anti-Mumps virus Non-structural protein V Polyclonal Antibody

Product Details

Reactivity: Mumps virus
Molecular Weight: Actual: 36 kDa.

Applications

Sample:
Lane 1: Recombinant Mumps virus Non-structural protein V
Verified Activity: Primary: Anti-Mumps virus Non-structural protein V (TMAB-09110) at 1/1000 dilution
Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution
Predicted band size: kDa
Observed band size: 36 kDa
Application: WB
Recommended WB: 1:500-2000

Properties

Shipping: Shipping with blue ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481