

## Anti-MYBBP1A Antibody (5G953)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 149 kDa. Actual: 149 kDa.
Clone:	5G953
Purification:	Protein A purified

### Applications

Verified Activity:	Western blot analysis of HEK293 cell lysate. Using MYBBP1A (TMAB-09134) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IHC-P,IHC-Fr,ICC/IF,IF,FCM
Recommended	WB: 1:1000-2000; IHC-P: 1:100-200; IHC-Fr: 1:100-200; ICC/IF: 1:50-200; IF: 1:50-200; FCM: 1:50-100

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Myb binding protein 1A
Antigen Species:	Human
Gene ID:	2617
Uniprot ID:	P41250

### Research Background

May activate or repress transcription via interactions with sequence specific DNA-binding proteins.Repression may be mediated at least in part by histone deacetylase activity.Acts as a corepressor and in concert with CRY1, represses the transcription of the core circadian clock component PER2.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481