

Anti-N myc interactor Antibody (80949)

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 35 kDa. Actual: 39 kDa.
Clone:	80949
Purification:	Protein A purified

Applications

Verified Activity:	25 µg total protein per lane of various lysates (see on figure) probed with N myc interactor monoclonal antibody, unconjugated (TMAB-09200) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IHC-P,IHC-Fr,ICC/IF,IF,IP
Recommended	WB: 1:500-2000; IHC-P: 1:100-200; IHC-Fr: 1:100-200; ICC/IF: 1:50-200; IF: 1:50-200; IP: 1:20-50

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human N myc interactor
Antigen Species:	Human
Gene ID:	9111
Uniprot ID:	Q13287

Research Background

NMYC interactor (NMI) encodes a protein that interacts with NMYC and CMYC (two members of the oncogene Myc family), and other transcription factors containing a Zip, HLH, or HLH Zip motif.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481