

## Anti-NCR3 Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Rat,Dog,Horse,Rabbit)
Molecular Weight:	Theoretical: 20 kDa.
Purification:	Protein A purified

### Applications

Verified Activity:	<p>1. Black line : Positive blank control (HL60); Negative blank control (HUVEC) Green line : Primary Antibody (Rabbit Anti-NCR3 antibody (TMAB-09290)) Orange line: Isotype Control Antibody (Rabbit IgG). Blue line : Secondary Antibody (Goat anti-rabbit IgG-AF647) HL60 (Positive) and Huvec (Negative control) Cells (black) were incubated in 5% BSA blocking buffer for 30 min at room temperature. Cells were then stained with NCR3 antibody (TMAB-09290) at 1:50 dilution in blocking buffer and incubated for 30 min at room temperature, washed twice with 2% BSA in PBS, followed by secondary antibody (blue) incubation for 40 min at room temperature. Acquisitions of 20,000 events were performed. Cells stained with primary antibody (green), and isotype control (orange).</p> <p>2. Blank control: HL60. Primary Antibody (green line): Rabbit Anti-NCR3 antibody (TMAB-09290) Dilution: 3 <math>\mu\text{g}</math> / <math>10^6</math> cells; Isotype Control Antibody (orange line): Rabbit IgG. Secondary Antibody: Goat anti-rabbit IgG-PE Dilution: 1 <math>\mu\text{g}</math> /test. Protocol The cells were incubated in 5% BSA to block non-specific protein-protein interactions for 30 min at room temperature. Cells stained with Primary Antibody for 30 min at room temperature. The secondary antibody used for 40 min at room temperature. Acquisition of 20,000 events was performed.</p>
Application:	FCM
Recommended	FCM: 2 $\mu\text{g}$ /Test

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen: KLH conjugated synthetic peptide: human CD337  
Antigen Species: Human  
Gene ID: 259197  
Uniprot ID: O14931

---

### Research Background

NCR3 is a cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis. It is a member of the Ig superfamily and may cooperate with NKp46 and NKp44 in the induction of cytotoxicity against a variety of target cells and is selectively expressed by all NK cells, both freshly isolated and cultured in IL2, thus representing an optimal marker for NK cell identification. Similar to NKp46 and CD16, NKp30 associated with CD3 zeta chains becomes tyrosine phosphorylated following cell treatment with pervanadate. Its extracellular portion is characterized by a single domain of the V-type and by a region rich in hydrophobic amino acids, potentially involved in protein/protein interactions that connects the Ig V-like domain with the transmembrane portion. It is selectively expressed by all resting and activated NK cells and weakly expressed in spleen.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481