

## Anti-NDRG3 Antibody (1Z937)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 41 kDa. Actual: 41 kDa.
Clone:	1Z937
Purification:	Protein A purified

### Applications

Verified Activity:	Western blot analysis of Raji cell lysate. Using NDRG3 (TMAB-09302) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IHC-P,IHC-Fr,IF
Recommended	WB=1:500-2000,IHC-P=1:100-200,IHC-F=1:100-200,IF=1:100-200

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human NDRG3
Antigen Species:	Human
Gene ID:	57446
Uniprot ID:	Q9UGV2

### Research Background

Predicted to be involved in signal transduction. Located in cytoplasm.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481