

## Anti-NDUFAF1 Antibody (612)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 38 kDa. Actual: 33 kDa.
Clone:	612
Purification:	Protein A purified

### Applications

Verified Activity:	Western blot analysis of 293T cell lysate. Using NDUFAF1 (TMAB-09310) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IHC-P,IHC-Fr,ICC/IF,IF,IP,FCM
Recommended	WB: 1:1000-2000; IHC-P: 1:100-200; IHC-Fr: 1:100-200; ICC/IF: 1:50-200; IF: 1:50-200; IP: 1:20-50; FCM: 1:50-100

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human NDUFAF1
Antigen Species:	Human
Gene ID:	51103
Uniprot ID:	Q9Y375

### Research Background

As part of the MCIA complex, involved in the assembly of the mitochondrial complex I.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481