

## Anti-NEDD4-2 Antibody (3Y211)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 112 kDa. Actual: 135 kDa.
Clone:	3Y211
Purification:	Protein A purified

### Applications

Verified Activity:	25 µg total protein per lane of various lysates (see on figure) probed with NEDD4-2 monoclonal antibody, unconjugated (TMAB-09336) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IP
Recommended	WB: 1:500-2000; IP: 1:20-50

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human NEDD4L
Antigen Species:	Human
Gene ID:	23327
Uniprot ID:	Q96PU5

### Research Background

E3 ubiquitin-protein ligase that mediates the polyubiquitination of lysine and cysteine residues on target proteins and is thereby implicated in the regulation of various signaling pathways including autophagy, innate immunity or DNA repair.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481