

Anti-Netrin 4 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Mouse (predicted:Human,Rat,Dog,Cow,Horse,Rabbit)
Molecular Weight:	Theoretical: 70 kDa. Actual: 63 kDa.
Purification:	Protein A purified

Applications

Sample:	Kidney (Mouse) Lysate at 40 µg Spleen (Mouse) Lysate at 40 µg Uterus (Mouse) Lysate at 40 µg
Verified Activity:	Heart (Mouse) Lysate at 40 µg Primary: Anti-Netrin 4 (TMAB-09366) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 70 kD Observed band size: 63 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human Netrin 4
Antigen Species:	Human
Gene ID:	59277
Uniprot ID:	Q9HB63

Research Background

May play an important role in neural, kidney and vascular development. Promotes neurite elongation from olfactory bulb explants.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481