

Anti-NFAM1 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat,Pig,Horse)
Molecular Weight:	Theoretical: 30 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	1. Blank control: U937. Primary Antibody (green line): Rabbit Anti-NFAM1 antibody (TMAB-09424) Dilution: 1 µg/Test; Secondary Antibody: Goat anti-rabbit IgG-AF488 Dilution: 0.5 µg/Test. Protocol The cells were incubated in 5% BSA to block non-specific protein-protein interactions for 30 min at room temperature. Cells stained with Primary Antibody for 30 min at room temperature. The secondary antibody used for 40 min at room temperature. Acquisition of 20,000 events was performed.
	2. Blank control: U937. Primary Antibody (green line): Rabbit Anti-NFAM1 antibody (TMAB-09424) Dilution: 1 µg/Test; Secondary Antibody: Goat anti-rabbit IgG-AF488 Dilution: 0.5 µg/Test. Protocol The cells were incubated in 5% BSA to block non-specific protein-protein interactions for 30 min at room temperature. Cells stained with Primary Antibody for 30 min at room temperature. The secondary antibody used for 40 min at room temperature. Acquisition of 20,000 events was performed.
Application:	FCM
Recommended	FCM: 1µg/Test

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: KLH conjugated synthetic peptide: human NFAM1
Antigen Species: Human
Gene ID: 150372
Uniprot ID: Q8NET5

Research Background

NFAM1 activates cytokine gene promoters such as the IL13 and TNF-alpha promoters. It contains an immunoreceptor tyrosine-based activation motif (ITAM) and is thought to regulate the signaling and development of B-cells.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481