

Anti-NOS2/iNOS Antibody (7E401)

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse)
Molecular Weight:	Theoretical: 131 kDa. Actual: 131 kDa.
Clone:	7E401
Purification:	Protein A purified

Applications

Verified Activity:	The HepG2 (treated with LPS (10 µg/mL, 6 h)) (H) cells were fixed with 4% PFA (10 min at r. T.) and then permeabilized with 90% ice-cold methanol for 20 min at -20°C, the cells then were incubated in 5% BSA to block non-specific protein-protein interactions (30 min at r. T.), followed by secondary antibody incubation for 40 min at room temperature. Primary Antibody (green): Rabbit Anti-iNOS antibody (TMAB-09566, 1:100); Isotype Control (orange): Rabbit IgG. Blank control (black): PBS. Acquisition of 20,000 events was performed.
Application:	FCM
Recommended	FCM: 1:50-100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human iNOS
Antigen Species:	Human
Gene ID:	4843
Uniprot ID:	P35228

Research Background

iNOS Produces nitric oxide (NO) which is a messenger molecule with diverse functions throughout the body. In macrophages, NO mediates tumoricidal and bactericidal actions.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481