

## Anti-OR5L1 Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Chimpanzee,Rhesus monkey,Gorilla)
Molecular Weight:	Theoretical: 36 kDa. Actual: 34 kDa.
Purification:	Protein A purified

## Applications

Verified Activity:	Sample: A431 Cell (Human) Lysate at 40 µg Primary: Anti-OR5L1 (TMAB-09793) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 36 kD Observed band size: 34 kD
Application:	WB
Recommended	WB: 1:500-2000

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human OR5L1
Antigen Species:	Human
Gene ID:	219437
Uniprot ID:	Q8NGL2

## Research Background

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481