

Anti-OTOL1 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat)
Molecular Weight:	Theoretical: 47 kDa. Actual: 49 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	1. Sample: LOVo (Human) Cell Lysate at 30 µg Primary: Anti-OTOL1 (TMAB-09831) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 47 kD Observed band size: 49 kD
	2. Sample: SW480 (Human) Cell Lysate at 30 µg Primary: Anti-OTOL1 (TMAB-09831) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 47 kD Observed band size: 49 kD
	Application: WB
	Recommended WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human OTOL1
Antigen Species:	Human
Gene ID:	131149
Uniprot ID:	A6NHN0

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481