

## Anti-Pantophysin Antibody (6C37)

## Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 29 kDa. Actual: 29-37 kDa.
Clone:	6C37
Purification:	Protein A purified

## Applications

Verified Activity:	<p>1. Paraformaldehyde-fixed, paraffin embedded Human Ovarian Cancer; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation Pantophysin Monoclonal Antibody, Unconjugated (TMAB-09968) at 1: 200 overnight at 4°C, followed by conjugation to the Goat Anti-Rabbit IgG H&amp;L Secondary Antibody-HRP and DAB staining.</p> <p>2. Paraformaldehyde-fixed, paraffin embedded Human Colon; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation Pantophysin Monoclonal Antibody, Unconjugated (TMAB-09968) at 1: 200 overnight at 4°C, followed by conjugation to the Goat Anti-Rabbit IgG H&amp;L Secondary Antibody-HRP and DAB staining.</p>
Application:	IHC-P,IHC-Fr,IF
Recommended	IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	A synthesized peptide: human SYPL1
Antigen Species:	Human
Gene ID:	6856
Uniprot ID:	Q16563

## Research Background

Predicted to be involved in chemical synaptic transmission. Located in extracellular exosome.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481