

Anti-PCBP1 Antibody (91764)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Molecular Weight:	Theoretical: 37 kDa. Actual: 37 kDa.
Clone:	91764
Purification:	Protein A purified

Applications

Verified Activity:	25 µg total protein per lane of various lysates (see on figure) probed with PCBP1 monoclonal antibody, unconjugated (TMAB-10037) at 1: 2000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IP
Recommended	WB: 1:500-2000; IP: 1:20-50

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human PCBP1
Antigen Species:	Human
Gene ID:	5093
Uniprot ID:	Q15365

Research Background

Single-stranded nucleic acid binding protein that binds preferentially to oligo dC. Together with PCBP2, required for erythropoiesis, possibly by regulating mRNA splicing (By similarity).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481