

## Anti-PER2 Polyclonal Antibody 2

### Product Details

Ig Type:	IgG
Reactivity:	Mouse,Rat (predicted:Human)
Molecular Weight:	Theoretical: 138 kDa.
Purification:	Protein A purified

### Applications

Verified Activity:	<p>1. Tissue/cell: rat skin tissue; 4% Paraformaldehyde-fixed and paraffin-embedded; Antigen retrieval: citrate buffer (0.01 M, pH 6.0), Boiling bathing for 15 min; Block endogenous peroxidase by 3% Hydrogen peroxide for 30 min; Blocking buffer (normal goat serum) at 37°C for 20 min; Incubation: Anti-PER2/Period circadian protein 2 Polyclonal Antibody, Unconjugated (TMAB-10178) 1: 200, overnight at 4°C, followed by conjugation to the secondary antibody and DAB staining</p> <p>2. Paraformaldehyde-fixed, paraffin embedded (Rat skeletal muscle); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (PER2) Polyclonal Antibody, Unconjugated (TMAB-10178) at 1:400 overnight at 4°C, followed by operating according to SP Kit (Rabbit) instructions and DAB staining.</p> <p>3. Paraformaldehyde-fixed, paraffin embedded (Mouse brain); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (PER2) Polyclonal Antibody, Unconjugated (TMAB-10178) at 1:400 overnight at 4°C, followed by operating according to SP Kit (Rabbit) instructions and DAB staining.</p>
Application:	IHC-P,IHC-Fr,IF
Recommended	IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen: KLH conjugated synthetic peptide: human PER2  
Antigen Species: Human  
Gene ID: 8864  
Uniprot ID: O15055

---

### Research Background

PER2, a mammalian homologue of the *Drosophila* period gene, shares a 40% homology with PER1 including the protein dimerization PAS domain. PER2 is a circadian regulator that may act as a transcription factor. It behaves as a negative element in circadian transcriptional loop. PER2 does not appear to bind DNA, suggesting indirect transcriptional inhibition. Expression oscillates with a 24 hour rhythm in the suprachiasmatic nucleus (SCN) and the whole eyes. Oscillations are maintained under constant darkness and are responsive to changes of the light/dark cycles. There is a 4 hour time delay between PER1 and PER2 oscillations. The expression rhythms appear to originate from retina.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481