

## Anti-PER-3 Antibody (7P192)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Molecular Weight:	Theoretical: 132 kDa. Actual: 130 kDa.
Clone:	7P192
Purification:	Protein A purified

### Applications

Verified Activity:	25 ug total protein per lane of various lysates (see on figure) probed with PER-3 monoclonal antibody, unconjugated (TMAB-10183) at 1: 1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB
Recommended	WB=1:500-2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human PER3
Antigen Species:	Human
Gene ID:	8863
Uniprot ID:	P56645

### Research Background

Corepressor targeting diverse transcription regulators. Functions in brown adipose tissue (BAT) differentiation.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481