

Anti-Perilipin A+B Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human, Mouse (predicted: Chicken, Dog, Sheep)
Molecular Weight:	Theoretical: 54-56 kDa. Actual: 50-56 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	1. Blank control: MCF7 (blue). Primary Antibody: Rabbit Anti-Perilipin A+B/AF647 Conjugated antibody (TMAB-10189-AF647), Dilution: 5µg in 100 µL 1X PBS containing 0.5% BSA; Isotype Control Antibody: Rabbit IgG/AF647 (orange), used under the same conditions.
	2. Blank control: U937 (fixed with 2% paraformaldehyde (10 min)). Primary Antibody: Rabbit Anti-Perilipin A+B antibody (TMAB-10189), Dilution: 0.2 µg in 100 µL 1X PBS containing 0.5% BSA; Isotype Control Antibody: Rabbit IgG (orange), used under the same conditions); Secondary Antibody: Goat anti-rabbit IgG-PE (white blue), Dilution: 1: 200 in 1 X PBS containing 0.5% BSA.
	3. This image was generously provided by James Fox, PhD. at the University of York. Immunostaining of perilipin-positive adipocytes (red) using Rabbit Anti-Perilipin A+B Polyclonal Antibody (TMAB-10189) in bone marrow of IL-7cre Rosa26-EYFP mice with DAPI (blue) nuclear counterstain.
	4. Sample: Stomach (Mouse) Lysate at 40 µg Primary: Anti-Perilipin A+B (TMAB-10189) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 54-56 kD Observed band size: 56 kD
	5. Sample: HepG2 Cell (Human) Lysate at 30 µg Primary: Anti-Perilipin A+B (TMAB-10189) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 54-56 kD Observed band size: 50 kD
Application:	WB, IHC-P, IHC-Fr, IF, FCM
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500; FCM: 0.2µg/Test

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: KLH conjugated synthetic peptide: human Perilipin A+B
Antigen Species: Human
Gene ID: 5346
Uniprot ID: O60240

Research Background

Perilipins, members of the PAT protein family (named after lipid droplet proteins Perilipin, Adipophilin, and TIP47) are found exclusively at the surface of lipid droplets in adipocytes and steroidogenic cells. They have been suggested to function as regulators of lipolysis and triacylglycerol storage within adipose tissue. Four distinct isoforms ranging from perilipin A (57 kDa) to perilipin D (26 kDa) have been identified and they share an identical amino terminal sequences, and contain 2-6 consensus protein kinase A (PKA) phosphorylation sites. Perilipin C and D have been detected only in steroidogenic cells. Perilipin A is the most abundant form on the lipid droplets of adipocytes. The phosphorylation of perilipin by PKA, which is accompanied by the phosphorylation and translocation of hormone-sensitive lipase from the cytosol to the lipid droplets, promotes lipolysis. There is evidence for the presence of perilipin A in atheroma plaques suggesting that the protein may be involved in the development oftherosclerosis by controlling as in adipocytes the hydrolysis of stored lipids.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481