

Anti-PRKAR1 Antibody (3T328)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Molecular Weight:	Theoretical: 43 kDa. Actual: 48 kDa.
Clone:	3T328
Purification:	Protein A purified

Applications

Verified Activity:	Western blot analysis of U-87MG cell lysate. Using PRKAR1 (TMAB-11738) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IHC-P,IHC-Fr,ICC/IF,IF
Recommended	WB: 1:1000-2000; IHC-P: 1:100-200; IHC-Fr: 1:100-200; ICC/IF: 1:50-200; IF: 1:50-200

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human PRKAR1A
Antigen Species:	Human
Gene ID:	5573
Uniprot ID:	P10644

Research Background

Regulatory subunit of the cAMP-dependent protein kinases involved in cAMP signaling in cells.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481