

Anti-PRL Antibody (2U760)

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 26 kDa. Actual: 22 kDa.
Clone:	2U760
Purification:	Protein A purified

Applications

Verified Activity:	1. Western blot analysis of Human fetal pituitary lysate. Using PRL (TMAB-11751) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min. 2. Paraformaldehyde-fixed, paraffin embedded Human pituitary; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with PRL Monoclonal Antibody, Unconjugated (TMAB-11751) at 1: 200 overnight at 4°C, followed by conjugation to the Goat Anti-Rabbit IgG H&L Secondary Antibody-HRP and DAB staining.
Application:	WB,IHC-P,IHC-Fr,IF
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Prolactin
Antigen Species:	Human
Gene ID:	5617
Uniprot ID:	P01236

Research Background

Prolactin acts primarily on the mammary gland by promoting lactation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481