

## Anti-PROX2 Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Mouse (predicted:Human,Rat,Dog,Pig,Cow,Horse,Rabbit,Sheep)
Molecular Weight:	Theoretical: 66 kDa. Actual: 66 kDa.
Purification:	Protein A purified

## Applications

Verified Activity:	Sample: Uterus (Mouse) Lysate at 40 µg Primary: Anti-PROX2 (TMAB-11814) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 66 kD Observed band size: 66 kD
Application:	WB
Recommended	WB: 1:500-2000

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human PROX2
Antigen Species:	Human
Gene ID:	283571
Uniprot ID:	Q3B8N5

## Research Background

The family of guanylin regulatory peptides, including guanylin and uroguanylin, are strongly expressed in intestinal mucosa and regulate intestinal fluid secretion during digestion. Guanylins are also involved in acid neutralization and the regulation of membrane-bound guanylate cyclase signaling molecules. Guanylin and uroguanylin are secreted primarily in the stomach, intestine, and colon. Uroguanylin is an endogenous activator of intestinal guanylate cyclase. It is a paracrine and/or autocrine regulator of intestinal water and salt transport. Uroguanylin stimulates intestinal guanylate cyclase through the same receptor binding region as the heat-stable enterotoxins. Uroguanylin is involved in the regulation of intestinal fluid and electrolyte transport.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481