

Anti-PRR16/DSC54 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human, Mouse (predicted: Rat, Chicken, Dog, Cow, Horse, Sheep)
Molecular Weight:	Theoretical: 33 kDa. Actual: 33 kDa.
Purification:	Protein A purified

Applications

Sample:	U-87Mg (Human) Lysate at 30 µg Cerebrum (Mouse) Lysate at 40 µg
Verified Activity:	Primary: Anti-DSC54 (TMAB-11827) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 33 kD Observed band size: 33 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human PRR16/DSC54
Antigen Species:	Human
Gene ID:	51334
Uniprot ID:	Q569H4

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481