

## Anti-PRSS2 Antibody (8M622)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Molecular Weight:	Theoretical: 26 kDa. Actual: 28 kDa.
Clone:	8M622
Purification:	Protein A purified

### Applications

Verified Activity:	Western blot analysis of Human pancreas lysate. Using PRSS2 (TMAB-11832) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB
Recommended	WB: 1:1000-2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Trypsin 2
Antigen Species:	Human
Gene ID:	5645
Uniprot ID:	P07478

### Research Background

Part of the insoluble cornified cell envelope (CE) of stratified squamous epithelia.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481