

Anti-PRSS42 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Mouse
Molecular Weight:	Theoretical: 32 kDa. Actual: 32 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Sample: Lung (Mouse) Lysate at 40 µg Primary: Anti-PRSS42 (TMAB-11836) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 32 kD Observed band size: 32 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: mouse PRSS42
Antigen Species:	Mouse
Gene ID:	235628
Uniprot ID:	Q8VIF2

Research Background

PRSS42 (Protease, Serine 42) is a Protein Coding gene. GO annotations related to this gene include serine-type endopeptidase activity. An important paralog of this gene is PRSS21.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481