

Anti-PTH1R Antibody (4V983)

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 66 kDa. Actual: 65 kDa.
Clone:	4V983
Purification:	Protein A purified

Applications

Verified Activity:	25 µg total protein per lane of various lysates (see on figure) probed with PTH1R monoclonal antibody, unconjugated (TMAB-11906) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IP
Recommended	WB: 1:500-2000; IP: 1:20-50

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human PTH1R
Antigen Species:	Human
Gene ID:	5745
Uniprot ID:	Q03431

Research Background

Receptor for parathyroid hormone and for parathyroid hormone-related peptide. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase and also a phosphatidylinositol-calcium second messenger system.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481