

Anti-PWWP2 Polyclonal Antibody

Product Details

| | |
|-------------------|---|
| Ig Type: | IgG |
| Reactivity: | Mouse,Rat (predicted:Human,Chicken,Horse,Zebrafish) |
| Molecular Weight: | Theoretical: 64 kDa. Actual: 64 kDa. |
| Purification: | Protein A purified |

Applications

| | |
|--------------------|--|
| Verified Activity: | 1. Sample: Testis (Mouse) Lysate at 40 µg Primary: Anti-PWWP2 (TMAB-11954) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 64 kD Observed band size: 63 kD |
| | 2. Sample: Fatal brain (Mouse) Lysate at 40 µg Primary: Anti-PWWP2 (TMAB-11954) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 64 kD Observed band size: 63 kD |
| Application: | WB |
| Recommended | WB: 1:500-2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|---|
| Immunogen: | KLH conjugated synthetic peptide: human PWWP2 |
| Antigen Species: | Human |
| Gene ID: | 170394 |
| Uniprot ID: | Q6NUJ5 |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481