

Anti-RAP80 Antibody (7N70)

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 80 kDa. Actual: 75 kDa.
Clone:	7N70
Purification:	Protein A purified

Applications

Verified Activity:	<p>1. 25 µg total protein per lane of various lysates (see on figure) probed with RAP80 monoclonal antibody, unconjugated (TMAB-12089) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.</p> <p>2. The MCF-7 (H) cells were fixed with 4% PFA (10 min at r. T.) and then permeabilized with 90% ice-cold methanol for 20 min at -20°C, the cells then were incubated in 5% BSA to block non-specific protein-protein interactions (30 min at r. T.), followed by secondary antibody incubation for 40 min at room temperature. Primary Antibody (green): Rabbit Anti-RAP80 antibody (TMAB-12089, 1:100); Isotype Control (orange): Rabbit IgG. Blank control (black): PBS. Acquisition of 20,000 events was performed.</p>
Application:	WB,IHC-P,IHC-Fr,ICC/IF,IF,FCM
Recommended	WB=1:500-2000,IHC-P=1:50-200,IHC-F=1:50-200,ICC/IF=1:50-200,IF=1:50-200,Flow-Cyt=1:50-100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human RAP80
Antigen Species:	Human
Gene ID:	51720
Uniprot ID:	Q96RL1

Research Background

Ubiquitin-binding protein. Specifically recognizes and binds 'Lys-63'-linked ubiquitin. Plays a central role in the BRCA1-A complex by specifically binding 'Lys-63'-linked ubiquitinated histones H2A and H2AX at DNA lesions sites, leading to target the BRCA1-BARD1 heterodimer to sites of DNA damage at double-strand breaks (DSBs).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481