

Anti-RBCK1 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Mouse (predicted:Human,Rat,Dog,Cow,Rabbit,Sheep)
Molecular Weight:	Theoretical: 56 kDa. Actual: 56 kDa.
Purification:	Protein A purified

Applications

Sample:	Uterus (mouse) Lysate at 40 µg Ovary (Mouse) Lysate at 40 µg
Verified Activity:	Primary: Anti-RBCK1/HBV associated factor 4 at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 56 kD Observed band size: 56 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human RBCK1/HBV associated factor 4
Antigen Species:	Human
Gene ID:	10616
Uniprot ID:	Q9BYM8

Research Background

The protein encoded by this gene is similar to mouse UIP28/UbcM4 interacting protein. Alternative splicing has been observed at this locus, resulting in distinct isoforms. [provided by RefSeq, Jul 2008]

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481