

## Anti-Renalase Antibody (2K659)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Molecular Weight:	Theoretical: 38 kDa. Actual: 43 kDa.
Clone:	2K659
Purification:	Protein A purified

### Applications

Verified Activity:	Western blot analysis of 293T cell lysate. Using Renalase (TMAB-12187) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IHC-P,IHC-Fr,IF
Recommended	WB: 1:1000-2000; IHC-P: 1:100-200; IHC-Fr: 1:100-200; IF: 1:100-200

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Renalase
Antigen Species:	Human
Gene ID:	55328
Uniprot ID:	Q5VYX0

### Research Background

Renalase is a flavin adenine dinucleotide-dependent amine oxidase that is secreted into the blood from the kidney.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481