

## Anti-RET Polyclonal Antibody 2

### Product Details

Ig Type:	IgG
Reactivity:	Mouse,Rat (predicted:Human,Dog,Rabbit)
Molecular Weight:	Theoretical: 122 kDa. Actual: 122 kDa.
Purification:	Protein A purified

### Applications

Verified Activity:	<p>1. Paraformaldehyde-fixed, paraffin embedded (Mouse brain); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (RET) Polyclonal Antibody, Unconjugated (TMAB-12200) at 1:400 overnight at 4° C, followed by operating according to SP Kit (Rabbit) instructions and DAB staining.</p> <p>2. Tissue/cell: rat brain tissue; 4% Paraformaldehyde-fixed and paraffin-embedded; Antigen retrieval: citrate buffer (0.01 M, pH 6.0), Boiling bathing for 15 min; Block endogenous peroxidase by 3% Hydrogen peroxide for 30 min; Blocking buffer (normal goat serum) at 37°C for 20 min;</p> <p>Incubation: Anti-RET Polyclonal Antibody, Unconjugated (TMAB-12200) 1: 200, overnight at 4°C, followed by conjugation to the secondary antibody and DAB staining</p> <p>3. Sample: Kidney (Mouse) Lysate at 40 µg</p> <p>Primary: Anti-RET (TMAB-12200) at 1/1000 dilution</p> <p>Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution</p> <p>Predicted band size: 122 kD</p> <p>Observed band size: 122 kD</p>
Application:	WB,IHC-P,IHC-Fr,IF
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human RET
Antigen Species:	Human
Gene ID:	5979
Uniprot ID:	P07949

### Research Background

This gene, a member of the cadherin superfamily, encodes one of the receptor tyrosine kinases, which are cell-surface molecules that transduce signals for cell growth and differentiation. This gene plays a crucial role in neural crest development, and it can undergo oncogenic activation in vivo and in vitro by cytogenetic rearrangement.

## A DRUG SCREENING EXPERT

---

Mutations in this gene are associated with the disorders multiple endocrine neoplasia, type IIA, multiple endocrine neoplasia, type IIB, Hirschsprung disease, and medullary thyroid carcinoma. Two transcript variants encoding different isoforms have been found for this gene. Additional transcript variants have been described but their biological validity has not been confirmed. [provided by RefSeq, Jul 2008]

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481