

Anti-RIN1 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Rat (predicted:Human,Mouse,Pig,Cow)
Molecular Weight:	Theoretical: 84 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Paraformaldehyde-fixed, paraffin embedded (rat brain); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (RIN1) Polyclonal Antibody, Unconjugated (TMAB-12258) at 1: 200 overnight at 4°C, followed by operating according to SP Kit (Rabbit) instructions and DAB staining.
Application:	IHC-P,IHC-Fr,IF
Recommended	IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human RIN1
Antigen Species:	Human
Gene ID:	9610
Uniprot ID:	Q13671

Research Background

RIN1 is a Ras effector protein which can affect Ras signaling in a variety of ways. It enhances signaling from ABL1 and ABL2 which are involved in regulation of cytoskeletal remodeling. It may function as a guanine nucleotide exchange factor in activation of Rab5a. Rin1 also competes with RAF1 for binding of the activated form of Ras. It may play a role inhibiting the modulation of neuronal plasticity in formation of aversive memory, and has been implicated as a tumor suppressor gene in breast cancer.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481