

Anti-RUNDC3B Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat,Dog,Pig,Horse,Rabbit)
Molecular Weight:	Theoretical: 53 kDa. Actual: 53 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Sample: DU145 (Human) Cell Lysate at 30 µg Primary: Anti-RUNDC3B (TMAB-12420) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 53 kD Observed band size: 53 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human RUNDC3B
Antigen Species:	Human
Gene ID:	154661
Uniprot ID:	Q96NLO

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481