

## Anti-SLC26A10 Polyclonal Antibody 2

### Product Details

Ig Type:	IgG
Reactivity:	Mouse (predicted:Human)
Molecular Weight:	Theoretical: 60 kDa. Actual: 60 kDa.
Purification:	Protein A purified

### Applications

Verified Activity:	Sample: Heart (Mouse) Lysate at 40 µg Primary: Anti-SLC26A10 (TMAB-12874) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 60 kD Observed band size: 60 kD
Application:	WB
Recommended	WB: 1:500-2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human SLC26A10
Antigen Species:	Human
Gene ID:	65012
Uniprot ID:	Q8NG04

### Research Background

SLC26A10 (Solute Carrier Family 26 Member 10) is a Protein Coding gene. GO annotations related to this gene include antiporter activity and sulfate transmembrane transporter activity. An important paralog of this gene is

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481