

Anti-SMYD2 Antibody (2V843)

Product Details

Ig Type:	IgG1, Kappa
Reactivity:	Human (predicted:Rat)
Clone:	2V843
Purification:	Protein G purified

Applications

Verified Activity:	Cell line: HepG2
	Fixative: 100% Ice-cold methanol
	Permeabilization: 0.1% Triton X-100
	Primary Ab dilution: 1:50
	Primary incubation condition: 4°C overnight
Application:	Secondary Ab: Goat Anti-Rabbit IgG
	Nuclear counter stain: DAPI (Blue)
Recommended	Comment: Color green is the positive signal for TMAB-12990
Recommended	WB,ICC/IF
Recommended	WB: 1:500-1000; ICC/IF: 1:100-500

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Gene ID:	56950
----------	-------

Research Background

SET domain-containing proteins, such as SMYD2, catalyze lysine methylation (Brown et al., 2006 [PubMed 16805913]).[supplied by OMIM, Nov 2008]

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481