

Anti-SPATA22 Polyclonal Antibody

Product Details

| | |
|-------------------|--|
| Ig Type: | IgG |
| Reactivity: | Human (predicted:Mouse,Pig,Cow,Horse,Rabbit) |
| Molecular Weight: | Theoretical: 41 kDa. Actual: 41 kDa. |
| Purification: | Protein A purified |

Applications

| | |
|--------------------|--|
| Verified Activity: | Sample: 293T (Human) Cell Lysate at 30 µg Primary: Anti-SPATA22 (TMAB-13085) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 41 kD Observed band size: 41 kD |
| Application: | WB |
| Recommended | WB: 1:500-2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|---|
| Immunogen: | KLH conjugated synthetic peptide: human SPATA22 |
| Antigen Species: | Human |
| Gene ID: | 84690 |
| Uniprot ID: | Q8NHS9 |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481