

## Anti-Syntenin 1 Antibody (6T29)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 32 kDa. Actual: 30 kDa.
Clone:	6T29
Purification:	Protein A purified

### Applications

Verified Activity:	25 ug total protein per lane of various lysates (see on figure) probed with Syntenin 1 monoclonal antibody, unconjugated (TMAB-13341) at 1: 1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	FCM,ICC/IF,IF,IHC-Fr,IHC-P,IP,WB
Recommended	FCM=1:50-100; ICC/IF=1:50-200; IF=1:50-200; IHC-Fr=1:100-200; IHC-P=1:100-200; IP=1:20-50; WB=1:500-2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Syntenin 1
Antigen Species:	Human
Gene ID:	6386
Uniprot ID:	O00560

### Research Background

Multifunctional adapter protein involved in diverse array of functions including trafficking of transmembrane proteins, neuro and immunomodulation, exosome biogenesis, and tumorigenesis.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481