

Anti-TMEM56/TLCD4 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse)
Molecular Weight:	Theoretical: 30 kDa. Actual: 30 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Western blot analysis of lysates from Jurkat cells, primary antibody was diluted at 1:1000, 4° over night
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human TLCD4
Antigen Species:	Human
Gene ID:	148534

Research Background

Predicted to be involved in lipid homeostasis. Predicted to be integral component of membrane. Predicted to be active in endoplasmic reticulum. [provided by Alliance of Genome Resources, Apr 2022]

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481