

## Anti-TNFRSF10A Polyclonal Antibody 3

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 48 kDa. Actual: 23 kDa.
Purification:	Protein A purified

### Applications

Verified Activity:	Sample: Lane 1: Recombinant human DR4 protein, C-His Primary: Anti-TNFRSF10A (TMAB-13633) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 48 kDa Observed band size: 23 kDa
Application:	WB
Recommended	WB: 1:500-2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: human DR4 protein
Antigen Species:	Human
Gene ID:	8797
Uniprot ID:	O00220

### Research Background

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is activated by tumor necrosis factor-related apoptosis inducing ligand (TNFSF10/TRAIL), and thus transduces cell death signal and induces cell apoptosis. Studies with FADD-deficient mice suggested that FADD, a death domain containing adaptor protein, is required for the apoptosis mediated by this protein. [provided by RefSeq, Jul 2008]

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481