

Anti-VSTM5 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Mouse (predicted:Human,Rat,Pig,Cow,Horse,Sheep)
Molecular Weight:	Theoretical: 19 kDa. Actual: 17 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Sample: stomach (Mouse) Lysate at 40 µg Primary: Anti-VSTM5 (TMAB-14146) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 19 kD Observed band size: 17 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human VSTM5
Antigen Species:	Human
Gene ID:	387804
Uniprot ID:	A8MXK1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481