

Anti-VSV-G tag Antibody (6Q246)

Product Details

Ig Type:	IgG
Reactivity:	Species independent
Molecular Weight:	Actual: 45 kDa.
Clone:	6Q246
Purification:	Protein G purified

Applications

Sample:	Lane1: VSV-G-Tagged Fusion Protein Overexpression E. Coli Lysate at 2 µg Lane2: VSV-G-Tagged Fusion Protein Overexpression E. Coli Lysate at 2 µg
Verified Activity:	Primary: Anti-VSV-G tag (TMAB-14147) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution Predicted band size: 41 kD Observed band size: 45 kD
Application:	ELISA,FCM,ICC/IF,IF,IHC-Fr,IHC-P,WB
Recommended	ELISA=1:1000-5000; FCM=1 µg/Test; ICC/IF=1:100-500; IF=1:200-1000; IHC-Fr=1:200-1000; IHC-P=1:200-1000; WB=1:1000-5000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated VSV-G tag
------------	--------------------------

Research Background

Vesicular stomatitis virus (VSV), an enveloped RNA virus from the Rhabdoviridae family, is released from the plasma membrane of host cells by a process called budding. The glycoprotein (VSV-g) contains a domain in its extracellular membrane proximal stem that appears to be needed for efficient VSV budding.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481