

Anti-VZV gE Polyclonal Antibody 2

Product Details

Ig Type:	IgG
Reactivity:	Varicella Zoster Virus
Molecular Weight:	Theoretical: 65 kDa. Actual: 55 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Sample: Lane 1: Recombinant Varicella-zoster virus gE protein, N-His Primary: Anti-VZV gE (TMAB-14152) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 65 kDa Observed band size: 55 kDa
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Varicella-zoster virus gE protein
------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481