

Anti-YAP1 Polyclonal Antibody 2

Product Details

Ig Type:	IgG
Reactivity:	Human, Mouse (predicted: Rat, Chicken, Pig, Rabbit, Zebrafish, Sheep)
Molecular Weight:	Theoretical: 55 kDa. Actual: 68 kDa.
Purification:	Protein A purified

Applications

Sample:	293T (Human) Cell Lysate at 30 µg Testis (Mouse) Lysate at 40 µg
Verified Activity:	Primary: Anti-YAP1 (TMAB-14226) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 65-70 kD Observed band size: 68 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human YAP1
Antigen Species:	Human
Gene ID:	10413
Uniprot ID:	P46937

Research Background

This gene encodes the human ortholog of chicken YAP protein which binds to the SH3 domain of the Yes proto-oncogene product. This protein contains a WW domain that is found in various structural, regulatory and signaling molecules in yeast, nematode, and mammals, and may be involved in protein-protein interaction. [provided by RefSeq].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481