

Anti-ZmCCD8 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Zea mays
Molecular Weight:	Theoretical: 63 kDa. Actual: 63 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Protein: Corn stamens (Zea mays) lysate at 30 µg; Primary: rabbit Anti-ZmCCD8 (TMAB-14325) at 1:300; Secondary: HRP conjugated Goat-Anti-rabbit IgG at 1: 5000; Predicted band size: 63 kD Observed band size: 63 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: ZmCCD8
Antigen Species:	Zea mays
Gene ID:	100502522

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481