

Anti-ZNF740 Polyclonal Antibody

Product Details

| | |
|-------------------|---|
| Ig Type: | IgG |
| Reactivity: | Human,Mouse (predicted: Cow,Rabbit,Sheep) |
| Molecular Weight: | Theoretical: 22 kDa. Actual: 22 kDa. |
| Purification: | Protein A purified |

Applications

| | |
|--------------------|--|
| Verified Activity: | <p>1. Tissue/cell: human cervical carcinoma; 4% Paraformaldehyde-fixed and paraffin-embedded; Antigen retrieval: citrate buffer (0.01 M, pH 6.0), Boiling bathing for 15 min; Block endogenous peroxidase by 3% Hydrogen peroxide for 30 min; Blocking buffer (normal goat serum) at 37°C for 20 min;</p> <p>Incubation: Anti-ZNF740 Polyclonal Antibody, Unconjugated (TMAB-14416) 1:500, overnight at 4°C, followed by conjugation to the secondary antibody and DAB staining</p> <p>2. Sample: Kidney (Mouse) Lysate at 40 µg Primary: Anti-ZNF740 (TMAB-14416) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 22 kD Observed band size: 22 kD</p> |
| Application: | WB,IHC-P,IHC-Fr,IF |
| Recommended | WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|--|
| Immunogen: | KLH conjugated synthetic peptide: Human ZNF740 |
| Antigen Species: | Human |
| Gene ID: | 283337 |
| Uniprot ID: | Q8NDX6 |

Research Background

May be involved in transcriptional regulation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481