

Anti-ZNF764 Polyclonal Antibody 2

Product Details

Ig Type:	IgG
Reactivity:	Mouse (predicted:Rat)
Molecular Weight:	Theoretical: 45 kDa. Actual: 45 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Sample: Placenta (Mouse) Lysate at 40 µg Primary: Anti-ZNF764 (TMAB-14419) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 45 kD Observed band size: 45 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: mouse ZNF764
Antigen Species:	Mouse
Gene ID:	233893

Research Background

ZNF764 (Zinc Finger Protein 764) is a Protein Coding gene. Diseases associated with ZNF764 include 16p11.2 microdeletion and thyroid hormone resistance. Among its related pathways are Gene Expression and Gene Expression. An important paralog of this gene is ZNF785.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481