

PI3K-IN-60

Chemical Properties

CAS No. : 2504036-13-7


Formula: C₂₂H₂₇N₅O₃S

Molecular Weight: 441.55

Store at low temperature

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PI3K-IN-60 (Compound 34) is a selective PI3K γ inhibitor for studying autoimmune diseases and cancers such as multiple sclerosis, breast cancer, and pancreatic ductal adenocarcinoma.
Targets(IC50)	PI3K

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2647 mL	11.3237 mL	22.6475 mL
5 mM	0.4529 mL	2.2647 mL	4.5295 mL
10 mM	0.2265 mL	1.1324 mL	2.2647 mL
50 mM	0.0453 mL	0.2265 mL	0.4529 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Xiaolin Hao. Phosphatidylinositol 3-kinase inhibitors. WO2020210379A1. 2020-10-15.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481