


Olivanic acid

Chemical Properties

CAS No. : 64761-66-6
 Formula: C₁₁H₁₄N₂O₄S
 Molecular Weight: 270.31
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Olivanic acid is a β -lactam antibiotic with activity against both Gram-positive and Gram-negative bacteria, and it inhibits β -lactamase.
Targets(IC50)	Antibacterial, Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6995 mL	18.4973 mL	36.9946 mL
5 mM	0.7399 mL	3.6995 mL	7.3989 mL
10 mM	0.3699 mL	1.8497 mL	3.6995 mL
50 mM	0.074 mL	0.3699 mL	0.7399 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481